

Likit Partikül Sayım Cihazı

Designed with functionality in mind, the Lighthouse LS-60 Liquid Sampler is an integrated syringe sampler and optical particle counting system with an easy-to-use Windows™ based software interface. Following a quality standard of reliability and dependability, Lighthouse Worldwide Solutions developed the LS-60 Liquid Sampler using the latest in laser optical particle counter technology. With a Sample Flow Rate of 60mL/min, the LS-60 provides unprecedented measurement speed and accuracy. The LS-60 uses a high accuracy syringe pump flow control system and a high sample volume sensor. It is a perfect tool for multiple batch sampling applications. The graphical user interface makes the LS-60 easy to configure and operate. Data is automatically recorded and easily available in graph and tabular form in real time. Reports can be preconfigured and printed at the touch of a key